

STEPS FOR THE GAME



Step 1
Power ON the device



Step 4
Click randomize to change the color of wells.



Step 7
If correct, well turns green
If wrong, well turns red



Step 2
Open Your App & fill in the details of patient & therapist



Step 5
With all set, Click on START to start the game



Step 8
For each well the time taken & the total error is updated in app



Step 3
Connect the app to the device via Bluetooth (CMC-COHMET)



Step 6
The wells glow one by one and voice command will start.



Step 9
Finally Upload the score to the database.



CMC

COHMET
let's rehabilitate

Designed, Developed & Made in INDIA



CMC

COHMET
let's rehabilitate

Designed, Developed & Made in INDIA



CMC

COHMET
let's rehabilitate

Designed, Developed & Made in INDIA



CMC

COHMET
let's rehabilitate

Designed, Developed & Made in INDIA



from
Bioengineering Dept.
CMC, Vellore, Tamil Nadu

CMC | **COHMET**
let's rehabilitate

