




CMC | **COHMET**
let's rehabilitate




cohmet-cmc.netlify.app



Know You Device

- 1. Power button
- 2. Bluetooth indicator
- 3. Power indicator
- 4. Charging indicator
- 5. Charging socket
- 6. Circuit holder
- 7. Portable peg holder
- 8. Wood Casing
- 9. Peg
- 10. Well Light
- 11. Well
- 12. Well tray

STEPS TO USE COHMET



Step 1

Power ON the device




Step 6

The wells glows one by one and voice command will start.



Step 2

Open Your App & fill in the details of patient & therapist




Step 7

If correct, well turns green
If wrong, well turns red




Step 8

For each well the time taken is updated in the app



Step 9

Once all well are done, the total error is updated



Step 10


Finally Upload the score to the database.



Connect the app to the device via Bluetooth (CMC-COHMET)




Click randomize to change the color of wells.



With all set, Click on START to start the game

USE CASES

- Cognitive Development
- Motor Skills Development
- Rehabilitation Assessment



DOWNLOAD APP

CMC | COHMET
let's rehabilitate

Proudly Made In INDIA



FEATURES

- App Controlled
- Portable
- Dynamic Patterns
- Vibrant PEG colors
- Time Tracking
- Error Tracking
- Rechargeable
- Data Protection
- Inbuild Peg Holder
- Cloud Data Storage
- Different Modes
- Modifiable Attachments
- 3 Language support
- 3 LED Indicators
- Color Detection
- Sound Indicator
- Compact Design
- Affordable

designed and developed by Dept Of BIOENGINEERING CMC, VELLORE



CHARGING AND USE



- Charge the device for **4 Hours for 8 Hours** of running time.
- Don't charge overnight.
- While charging make sure RED light is on.
- Avoid Playing while charging.

Let's rehabilitate